AT1315 Gamma Beta Spectrometer



Two channel scintillation gamma and beta spectrometer is designed for simultaneous and selective detection of the following:



- ¹³⁷Cs, ¹³⁴Cs, ¹³¹I and ⁹⁰Sr radionuclide specific activity in natural samples
- Specific effective activity of ⁴⁰K, ²²⁶Ra, ²³²Th in construction materials

Can be used for rapid radioactive purity determination of standardized sample heat of metal.

Operating principle

Operating principle of AT1315 Gamma and beta spectrometer is based on measurement and transformation of gamma and beta radiation, detected by standalone detection units, into amplitude distributions, which are further transformed into digital code and saved in the detection unit memory.

The spectrometer can be delivered w/o beta channel according to customer request.



Gamma spectrometer

Applications

- Spectrometric and radiometric monitoring of radionuclide content in water, foods, agricultural raw materials, industrial, construction and forestry materials, environmental objects (soil, vegetation, etc.), metallurgical industry produce and products of other industries.
- Activity measurement of ¹³⁷Cs, ¹³⁴Cs, ¹³¹I,
 Sr, ⁴⁰K, ²²⁶Ra, ²³²Th, etc.

Features

- 1024-channel analogue-to-digital converter is integrated into smart probes
- Continuous automatic LED stabilisation of spectrometric path in measurement mode
- Calibration integrity monitoring by means of integrated radioisotope check sample with less than the minimum significant activity
- Computer spectra processing with maximum likelihood method
- Automatic recording of sample density
- Spectra metering with on-line visualisation
- Simultaneous spectra metering and processing
- Methodological support of measurements



AT1315 Gamma Beta Spectrometer

Specification

Detectors

Gamma channel Scintillator, NaI(TI) ø63x63 mm Beta channel Plastic scintillator, ø128x9 mm

Energy range

50 keV...3 MeV Gamma radiation Beta radiation 150 keV...3.5 MeV

Volumetric (specific) activity measuring

range without sample concentration (in Spectrometric and Radiometric modes)

¹³⁷Cs ^{40}K ²²⁶Ra

²³²Th 90Sr (In Radiometric mode only)

¹³¹I (In Spectrometric mode only) ¹³⁴Cs (In Spectrometric mode only)

Measured sample density range

Intrinsic relative error of activity

measurement with P = 0.95

Lower limit of 90 Sr measurement range

with sample concentration in conversion to "wet" sample

For potable water For milk, baby food For potatoes, corn, grain, agricultural raw materials 0.1 Bq/I 0.8 Bq/I

1...10⁶ Bq/I (Bq/kg)

3...10⁴ Bq/l (Bq/kg)

3...104 Bq/I (Bq/kg)

10...106 Bq/l (Bq/kg)

10...10⁵ Bq/l (Bq/kg)

6...105 Bq/l (Bq/kg)

±20% max.

0.2...1.6 g/cm3

20...2·10⁴ Bq/l (Bq/kg)

1.0 Bq/kg

<1%

5·104 s-1

<5%

1024

≥24 h

<15 min

≤75%

10°C...35°C

PC USB port

ø98x330 mm, 2 kg

ø138x323 mm, 2.5 kg

ø474x910 mm, 194 kg

Integral nonlinearity

Typical resolution at 662 keV (137Cs) 7.5%

Maximum input statistical load

Calibration scale instability <2%

during continuous service

Measurement instability during continuous service

Number of ADC channels

Continuous run time

Operation mode setup time

Working temperature range

Relative humidity with air temperature

Power supply

Overall dimensions, weight

≤30°C without condensation

Gamma detection unit

Beta detection unit Protection unit (with Beta detection unit)

Measurement vessels volume For "wet" samples

For concentrated samples


Marinelli beaker, 1 l; Flat vessel, 0.5 and 0.1 l Flat vessel, 0.2 I and 0.03 I

Design and specifications are subject to change without notice


Background measurement



Spectrometric mode



Radiometric mode



AT1315 Gamma beta spectrometer meets Safety standard requirements: IEC 61010-1:2001 EMC requirements: EN 55011:2009 IEC 61000-4-2:2008 IEC 61000-4-3:2008

Gamma beta spectrometer has the pattern approval certificates of Republic of Belarus, Russian Federation and Kazakhstan.





